

Bone Marrow Biopsy

Order Name: **None available**
Revision Date: 03/02/2016

TEST NAME		METHODOLOGY		LOINC CODE
Bone Marrow Biopsy				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	<p>Instructions:</p> <p>1. Prior to collection notify laboratory (405-307-1100) to schedule a lab tech to be present during procedure.</p> <p>2. Bone marrow specimens are collected by the physician and assisted by lab tech.</p> <p>3. The labeled specimen and slides are brought to Histology unfixed and in a syringe (if during open hours) accompanied by a completed tissue slip.</p> <p>4. If after hours, the specimen is placed in 10% neutral buffered formalin.</p> <p>Note: A pathology/diagnostic report and a brief history are required for processing.</p> <p>Storage temperature: Ambient</p> <p>Do not Refrigerate or Freeze specimen</p>			
Reference Range	An interpretive report will be provided.			
Methodology	Includes tissue processing, hematoxylin-and-eosin stain, iron stain, and additional stains if necessary for interpretation. A pathology consultation is also included and is charged separately.			

GENERAL INFORMATION	
Testing Schedule	Monday through Friday
CPT Code(s)	88305-Surgical pathology, gross and microscopic 88311-Decalcification
Lab Section	NRLS-Histology