

Mycoplasma pneumoniae Antibody, IgM

Order Name: **MYCO**
Revision Date: 01/12/2013

TEST NAME		METHODOLOGY		LOINC CODE
Mycoplasma pneumoniae Antibody, IgM				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Draw blood in a serum gel tube(s). Spin down and send 1 mL (minimum volume: 0.2 mL) of serum refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of draw, collector initials, and test(s) being ordered.			
Reference Range	Negative			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday TAT: 2 hours
CPT Code(s)	86738
Lab Section	NRLS-Core Serology