

Fusarium moniliforme Allergen

Order Name: **FUSARIUM**

Test Number: 5616830

Revision Date: 02/11/2013

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------------------|-------------|------------|
| Fusarium moniliforme Allergen | ImmunoCAP | 15722-2 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 0.4 mL (0.2 mL) | Serum | Clot Activator SST | Room Temperature |
| Instructions | Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3) | | | |

| GENERAL INFORMATION | |
|---------------------|---------------|
| Expected TAT | 3-5 days |
| Notes | |
| CPT Code(s) | 86003 |
| Lab Section | Reference Lab |