

## Meconium Drug Screen (5 Panel)

Order Name: **DRUG F SCR**

Test Number: 4066250

Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Meconium Drug Screen (5 Panel)	Immunoassay (IA)	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5g (1g)	See Instructions	See Instructions	Refrigerated
Instructions	<p><b>CHAIN OF CUSTODY TEST !</b>  <b>Specimen Type:</b> Meconium Collection Kit (PeopleSoft#96663) OR clean, leakproof polypropylene container.  <b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.  <b>Specimen Collection:</b> Specimens from different voidings may be pooled if necessary  <b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#4034</p>			

### GENERAL INFORMATION

Expected TAT	3-6 Days
Clinical Use	<p>All fecal drug screens are analyzed by immunoassay. Each positive result is confirmed by GC/MS or LC/MS/MS to the limit of detection.</p> <p><a href="#">Meconium Drug Screen Confirmations</a></p>
Performing Labcorp Test Code	809373
Notes	Labcorp Test Code: 809373
CPT Code(s)	80307
Lab Section	Reference Lab