

## Follicular Lymphoma, IGH/BCL2, t(14;18) by FISH

Order Name: **FOLLIC FSH**

Test Number: 9116975

Revision Date: 10/01/2022

| TEST NAME                                       | METHODOLOGY                        | LOINC CODE |
|---|------------------------------------|------------|
| Follicular Lymphoma, IGH/BCL2, t(14;18) by FISH | Fluorescence in Situ Hybridization |            |

### SPECIMEN REQUIREMENTS

| Specimen     | Specimen Volume (mL)  | Specimen Type | Specimen Container                 | Transport Environment |
|--------------|---|---------------|------------------------------------|-----------------------|
| Preferred    | 5 mL (1 mL)   | Bone Marrow   | Sodium Heparin (Green Top, No-Gel) | Room Temperature      |
| Alternate 1  | 5 mL (3 mL)   | Whole Blood   | Sodium Heparin (Green Top, No-Gel) | Room Temperature      |
| Instructions | Send specimen ASAP, Keep at room temperature! (DO NOT FREEZE). Frozen samples will be rejected. |               |                                    |                       |

### GENERAL INFORMATION

|                  |   |
|------------------|---|
| Testing Schedule | Mon-Fri   |
| Expected TAT     | 3-5 Days  |
| Clinical Use     | Useful to diagnose follicular lymphoma (> 75%) and some diffuse large B-cell lymphomas. |
| CPT Code(s)      | 88377   |
| Lab Section      | Reference Lab   |