

Celiac Testing at Prometheus

Order Name: **PRO CELIAC**

Test Number: 3613275

Revision Date: 12/08/2016

| TEST NAME                    | METHODOLOGY | LOINC CODE |
|------------------------------|-------------|------------|
| Celiac Testing at Prometheus |             |            |

| SPECIMEN REQUIREMENTS |                                                                                                                                                                                                                                                                                                                                                                          |                            |                                            |                       |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|--------------------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)                                                                                                                                                                                                                                                                                                                                                    | Specimen Type              | Specimen Container                         | Transport Environment |
| Preferred             | See Instructions                                                                                                                                                                                                                                                                                                                                                         | EDTA Whole Blood and Serum | EDTA (lavender top) and Clot Activator SST | Refrigerated          |
| Instructions          | Collect both serum and whole blood. 2mL serum (red/grey or red no-gel) and 5mL whole blood (EDTA).<br>Indicate desired test code: <ul style="list-style-type: none"><li>6360 - Celiac Plus</li><li>6260 - Celiac Genetics HLA DQ2/DQB</li><li>1155 - Celiac Serology</li></ul> Specimen Stability: Room temperature: 7 days, Refrigerated: 30 days, Frozen: Unacceptable |                            |                                            |                       |

| GENERAL INFORMATION |                           |
|---------------------|---------------------------|
| Expected TAT        | 3 days                    |
| Notes               | Reference Lab: Prometheus |
| Lab Section         | Reference Lab             |