

Helicobacter pylori Antibody IgG

Order Name: **H PYLORI Ab IgG**

Test Number: 5577751

Revision Date: 10/30/2019

| TEST NAME | METHODOLOGY | LOINC CODE |
|----------------------------|-------------------------|------------|
| H Pylori IgG | Chemiluminescence Assay | 5176-3 |
| H Pylor IgG Interpretation | | 16126-5 |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
|--------------|---|---------------|--------------------|-----------------------|
| Preferred | 1 mL (0.5mL) | Serum | Clot Activator SST | Refrigerated |
| Instructions | Testing should only be performed on adults (18 or older) Specimen Stability: Room temperature: Unacceptable, Refrigerated: 7 days, Frozen: 1 year | | | |

GENERAL INFORMATION

| | |
|------------------|--|
| Testing Schedule | Mon - Sat |
| Expected TAT | 1-2 Days |
| Clinical Use | Serology for Helicobacter pylori may be a useful noninvasive screening test for H. pylori infection. |
| CPT Code(s) | 86677 |
| Lab Section | Immunology - Serology |