

Deoxy pyridinoline Crosslinks (Pyrilinks-D)

Order Name: **Pyrilinks-D**

Test Number: 6904991

Revision Date: 01/28/2025

TEST NAME	METHODOLOGY	LOINC CODE
Deoxy pyridinoline Crosslinks (Pyrilinks-D)	Immunochemiluminometric (ICMA)	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, First Void	Sterile Urine container	Refrigerated

Instructions

Specimen Type: Plastic urine container
Specimen Storage: Refrigerated
Specimen Collection: Collection before 10 AM is recommended (eg, first morning void). Do NOT add preservative.
Special Instructions: The account must submit the patient's Social Security number to monitor.
Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days

GENERAL INFORMATION

Expected TAT	4-5 days
Performing Labcorp Test Code	511105
CPT Code(s)	82570, 82523
Lab Section	Reference Lab