

Mirtazapine (Remeron), Serum or Plasma

Order Name: **Mirtazapine**

Test Number: 6905779

Revision Date: 10/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Mirtazapine (Remeron), Serum or Plasma	High-pressure liquid chromatography (HPLC) with fluorometric detection	17283-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Instructions	Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube. Special Instructions: Testing referred to MedTox Laboratories Inc MX#434			

GENERAL INFORMATION	
Expected TAT	10-14 days
CPT Code(s)	80299
Lab Section	Reference Lab