

## Donath-Landsteiner Test

Order Name: **Donath-Land**

Test Number: 7306200

Revision Date: 09/24/2019

TEST NAME	METHODOLOGY	LOINC CODE
Donath-Landsteiner Test	Hemolysis	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (7 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Blood should be allowed to clot at 37°C, and the serum separated from the clot at this temperature. Once the serum is separated from the cells, the serum can be stored and transported refrigerated.			

### GENERAL INFORMATION

Testing Schedule	Mon-Fri (7:00am - 3:30pm)
Expected TAT	1 Day
Clinical Use	Used in the diagnosis of Paroxysmal Cold Hemoglobinuria (PCH).
CPT Code(s)	86940 86941
Lab Section	Blood Bank