

## Pain Panel, General Urine Drug Screen

Order Name: **PAIN PNLG**

Test Number: 4000314

Revision Date: 09/17/2018

TEST NAME	METHODOLOGY	LOINC CODE
Pain Panel, General Urine Drug Screen	Liquid Chromatography/Mass Spectroscopy	
Adulterants Screen Urine		
Prescribed Medication List - Prompt		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Container or Monovette	Refrigerated
Instructions	Keep refrigerated. Specimen Stability: Ambient 7 days. Refrigerated 30 days. Frozen 6 months Please provide accurate prescription list to lab.			

Labcorp Oklahoma, Inc.  
Test Directory

**GENERAL INFORMATION**

**Testing Schedule** Mon - Fri

**Expected TAT** 3-5 days

**Clinical Use** Used to evaluate medication compliance in pain management.  
Secondary testing is generally not indicated but is available if confirmation or quantitation is desired.

Screen (assay):	Confirmation (Order):
Acetaminophen	(4000126) Acetaminophen Confirmation
Amphetamines	(4000415) Amphetamine Class Confirmation
Barbiturates	(4000455) Barbiturate Class Confirmation
Benzodiazepines	(4000425) Benzodiazepine Class Confirmation
Buprenorphine	(4000127) Buprenorphine and Opiate Analog Confirmation
Carisoprodol	(4000130) Carisoprodol Confirmation
Fentanyl	(4000129) Fentanyl and Analog Confirmation
Gabapentin	(4000132) Gabapentin Confirmation
Methadone	(4000131) Methadone Confirmation
Opiates	(4000108) Opiate Class Confirmation
Pregabalin	(4000136) Pregabalin Confirmation
Propoxyphene	(4000133) Propoxyphene Confirmation
Tapentadol	(4000135) Tapentadol Confirmation
Tramadol	(4000137) Tramadol Confirmation
Tricyclics	(4000128) Tricyclic Confirmation
Zolpidem	(4000144) Zolpidem Confirmation

**Notes** **Adulterants Screen:**  
pH, Urine  
Creatinine, Urine

**CPT Code(s)** 80307

**Lab Section** Chemistry LC-MS