

## MySwab Vaginosis Profile Swab - TMA

Order Name: **MySwab Vag**  
Test Number: 6987003  
Revision Date: 04/06/2020

TEST NAME	METHODOLOGY	LOINC CODE
Bacterial Vaginosis - Swab	Transcription-Mediated Amplification	46999-9
Candida Species - Swab	Transcription-Mediated Amplification	47000-5
Candida Glabrata - Swab	Transcription-Mediated Amplification	n/a

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>Swab</b>	<b>Swab</b>	<b>APTIMA Multitest Swab kit</b>	<b>Refrigerated</b>
<b>Instructions</b>	Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87801, 87481
Lab Section	Molecular Diagnostic