

MySwab Vaginosis Profile + Trichomonas - Swab - TMA

Order Name: **MySwab Vag TV**

Test Number: 6987004

Revision Date: 04/06/2020

| TEST NAME                  |  | METHODOLOGY                          | LOINC CODE |
|----------------------------|--|--------------------------------------|------------|
| Bacterial Vaginosis - Swab |  | Transcription-Mediated Amplification | 46999-9    |
| Candida Species - Swab     |  | Transcription-Mediated Amplification | 47000-5    |
| Candida Glabrata - Swab    |  | Transcription-Mediated Amplification | n/a        |
| Trichomonas Swab (Vaginal) |  | Transcription-Mediated Amplification | 70165-6    |

| SPECIMEN REQUIREMENTS |                                                                                                                                                                                                                           |               |                           |                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|---------------------------|-----------------------|
| Specimen              | Specimen Volume (min)                                                                                                                                                                                                     | Specimen Type | Specimen Container        | Transport Environment |
| Preferred             | Swab                                                                                                                                                                                                                      | Swab          | APTIMA Multitest Swab kit | Refrigerated          |
| Instructions          | Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube.<br>Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A |               |                           |                       |

| GENERAL INFORMATION |                       |
|---------------------|-----------------------|
| Testing Schedule    | Mon-Fri               |
| Expected TAT        | 1-2 days after set-up |
| CPT Code(s)         | 87801, 87481, 87661   |
| Lab Section         | Molecular Diagnostic  |