

Mumps Virus PCR - Swab

Order Name: **Mumps PCR**

Test Number: 6905087

Revision Date: 09/26/2024

| TEST NAME              | METHODOLOGY               | LOINC CODE |
|------------------------|---------------------------|------------|
| Mumps Virus PCR - Swab | Polymerase Chain Reaction | 47532-7    |

| SPECIMEN REQUIREMENTS |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |               |                             |                       |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------------------|-----------------------|
| Specimen              | Specimen Volume (min)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | Specimen Type | Specimen Container          | Transport Environment |
| Preferred             | 1 swab                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Buccal swab   | Viral Transport Media (VTM) | Refrigerated          |
| Instructions          | <p><b>Specimen:</b> Buccal swab</p> <p><b>Container:</b> Viral Transport Media (VTM) or Universal Transport Media (UTM)</p> <p><b>Collection:</b> Collect a buccal swab specimen by massaging the parotid gland for 30 seconds. Using a synthetic swab, collect the sample by sweeping the swab between the cheek and gum from the upper molar to the lower molar for 10-15 seconds, while turning the swab. Sweep both sides with a single swab and place the swab into the VTM or UTM tube. Ensure the VTM or UTM tube is capped tightly to avoid leakage.</p> <p><b>Stability Requirements:</b> Ambient N/A; Refrigerated 7 days; Frozen 14 days</p> <p><b>Cause for Rejection:</b> Improperly submitted specimens: Specimens not received in an approved VTM/UTM; swabs received in Liquid Amies (E-Swabs) or other bacterial swabs; specimens received without a swab; grossly leaking specimens; specimens received outside the stated stability.</p> |               |                             |                       |

| GENERAL INFORMATION          |                     |
|------------------------------|---------------------|
| Expected TAT                 | 3-4 Days from setup |
| Performing Labcorp Test Code | 140283              |
| CPT Code(s)                  | 87798               |
| Lab Section                  | Reference Lab       |