

Procollagen I N-Terminal Propeptide

Order Name: Procollagen NTP

Test Number: 6906965

Revision Date: 12/12/2022

| TEST NAME                           | METHODOLOGY                   | LOINC CODE |
|-------------------------------------|-------------------------------|------------|
| Procollagen I N-Terminal Propeptide | Quantitative Radioimmunoassay | 47255-5    |

| SPECIMEN REQUIREMENTS |                                                                                                                                                                                                                                                                                                                                                                                                                               |               |                                  |                       |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|----------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)                                                                                                                                                                                                                                                                                                                                                                                                         | Specimen Type | Specimen Container               | Transport Environment |
| Preferred             | 1 mL (0.4 mL)                                                                                                                                                                                                                                                                                                                                                                                                                 | Serum         | Clot Activator SST               | Refrigerated          |
| Alternate 1           | 1 mL (0.4 mL)                                                                                                                                                                                                                                                                                                                                                                                                                 | Serum         | Clot Activator (Red Top, No-Gel) | Refrigerated          |
| Instructions          | <p><b>Notes:</b> 0.4 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> Red-top tube or gel-barrier tube</p> <p><b>Specimen Storage:</b> Refrigerated</p> <p><b>Specimen Collection:</b> If red-top tube is used, transfer separated serum to a plastic transport tube.</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p> |               |                                  |                       |

| GENERAL INFORMATION          |                           |
|------------------------------|---------------------------|
| Expected TAT                 | 4 - 6 days                |
| Performing Labcorp Test Code | 140850                    |
| Notes                        | Labcorp Test Code: 140850 |
| CPT Code(s)                  | 83519                     |
| Lab Section                  | Reference Lab             |