

## Peripheral Blood for Mycosis Fungoides by Flow

Order Name: **Mycosis Flow**  
Test Number: 5582601  
Revision Date: 07/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Peripheral Blood for Mycosis Fungoides by Flow	Flow cytometry	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	4.5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Instructions	<b>Deliver to laboratory (flow cytometry) ASAP.</b> DO NOT Centrifuge or Refrigerate. Collect Monday through Friday only! <b>Specimen must be received in the lab section by Friday afternoon or collect sample on Monday.</b> Specimen stability: 48hrs Room Temperature.			

### GENERAL INFORMATION

Testing Schedule	Tue-Sat
Expected TAT	2 days
Notes	Results from this test will be reported in a pathology report: <b>Order Code: [7120723] Flow Cytometry Final Report</b> <b>Result Code: [7120723] Flow Cytometry Final Report</b>
CPT Code(s)	88184, 88185X18, 88187
Lab Section	Flow Cytometry