

## Aspergillus (Galactomannan) Antigen

Order Name: **ASPAG**  
Revision Date: 10/08/2019

| TEST NAME                           | METHODOLOGY | LOINC CODE |
|-------------------------------------|-------------|------------|
| Aspergillus (Galactomannan) Antigen |             |            |

### SPECIMEN REQUIREMENTS

| Specimen               | Specimen Volume (min)   | Specimen Type | Specimen Container | Transport Environment |
|------------------------|---|---------------|--------------------|-----------------------|
| Preferred              | 1.5 mL (1 mL)   | Serum         | Serum Gel Tube     | Refrigerated          |
| <b>Instructions</b>    | Collection Instructions:  |               |                    |                       |
|                        | 1. Avoid exposure of specimen to atmosphere to prevent sample contamination from environment. |               |                    |                       |
|                        | 2. Centrifuge and send specimen in collection SST tube. <b>Do not aliquot or open tube.</b>   |               |                    |                       |
| <b>Reference Range</b> | <0.5 index  |               |                    |                       |
|                        | Reference values apply to all ages.   |               |                    |                       |
| <b>Methodology</b>     | Enzyme Immunoassay (EIA)  |               |                    |                       |

### GENERAL INFORMATION

|                         |  |
|-------------------------|--|
| <b>Testing Schedule</b> | Monday-Friday and Sunday   |
| <b>Expected TAT</b>     | 2 Days   |
| <b>Notes</b>            | Stability:<br>Refrigerated(preferred)- 14 Days<br>Frozen - 14 Days |
| <b>CPT Code(s)</b>      | 87305  |
| <b>Lab Section</b>      | NRLS-Mayo Medical Laboratories                                     |