

Estrone

Order Name: **E1**
Revision Date: 10/10/2019

TEST NAME		METHODOLOGY		LOINC CODE
Estrone				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.2 mL (0.8 mL)	Serum	Red Top	Room Temperature
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	See interpretive report.			
Methodology	Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)			

GENERAL INFORMATION	
Testing Schedule	Monday-Saturday; 1300
Expected TAT	3 Days
Notes	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.
CPT Code(s)	82679
Lab Section	NRLS-Mayo Medical Laboratories