

Alkaline Phosphatase Isoenzyme

Order Name: **ALKPHOSISO**
Revision Date: 11/08/2019

TEST NAME		METHODOLOGY		LOINC CODE
Alkaline Phosphatase Isoenzyme				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.5 mL)	Serum	Serum Gel Tube	Frozen
Alternate 1	1.0 mL (0.5 mL)	Serum	Red Top	Frozen
Instructions	Specimen must be divided into 2 tubes each containing 0.5 mL			
	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	ALP: Photometric, p-Nitrophenol Phosphate ALKE: Electrophoresis			

GENERAL INFORMATION	
Expected TAT	3 Days
Notes	Specimen Stability: Frozen (preferred) - 60 Days Ambient - 7 Days Refrigerated - 7 Days
CPT Code(s)	Alkaline Phosphatase, Serum 84075 Alkaline Phosphatase Isoenzymes, Serum 84080
Lab Section	NRLS-Mayo Medical Laboratories