

HLA-B27

Order Name: **LY27B**
Revision Date: 10/21/2019

TEST NAME		METHODOLOGY		LOINC CODE
HLA-B27				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	6 mL (1 mL)	Whole Blood	Lavender Top K2 EDTA	Room Temperature
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Flow Cytometry			

GENERAL INFORMATION	
Testing Schedule	Monday through Friday; 7:30-5 p.m. ; Saturday; 10 a.m.-6 p.m.
Notes	Specimen must arrive withing 96 hours of collection.
	Specimen Stability: Ambient (preferred) - 4 Days
CPT Code(s)	86812
Lab Section	NRLS-Mayo Medical Laboratories