

HLA-B27

Order Name: **LY27B**
Revision Date: 10/21/2019

TEST NAME	METHODOLOGY	LOINC CODE
HLA-B27		

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	6 mL (1 mL)	Whole Blood	Lavender Top K2 EDTA	Room Temperature
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Flow Cytometry			

GENERAL INFORMATION

Testing Schedule	Monday through Friday; 7:30-5 p.m. ; Saturday; 10 a.m.-6 p.m.
Notes	Specimen must arrive withing 96 hours of collection. Specimen Stability: Ambient (preferred) - 4 Days
CPT Code(s)	86812
Lab Section	NRLS-Mayo Medical Laboratories