

25-Hydroxyvitamin D2 and D3, Serum

Order Name: **25HDN**
Revision Date: 11/26/2019

TEST NAME	METHODOLOGY	LOINC CODE
25-Hydroxyvitamin D2 and D3, Serum		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.25 mL)	Serum	Red Top	Refrigerated
Alternate 1	0.5 mL (0.25 mL)	Serum	Serum Gel Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)			

GENERAL INFORMATION	
Testing Schedule	Monday through Friday continuous until 2:30 p.m.
Expected TAT	3-5 Days
Notes	Specimen Stability: Refrigerated (preferred)- 14 Days Frozen - 30 Days Ambient - 7 Days
CPT Code(s)	82306
Lab Section	NRLS-Mayo Medical Laboratories