## Labcorp Oklahoma, Inc. Test Directory

## **Bile Acids, Fractionated and Total**

Methodology

Order Name: BAFS

Revision Date: 12/02/2019

TEST NAME		M	ETHODOLOGY	LOINC CODE
Bile Acids, Fractionate	ed and Total			
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.3 mL)	Serum	Red Top	Refrigerated
Alternate 1	0.5 mL (0.3 mL)	Serum	Serum Gel Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	Analyte	Normal (nmol/mL)		
	Total Cholic Acid	< or =5.00		
	Total Chenodeoxycholic Acid	< or =6.00		
	Total Deoxycholic Acid	< or =6.00		
	Total Ursodeoxycholic Acid	< or =2.00		
	Total Bile Acids	< or =19.00		

GENERAL INFORMATION		
Testing Schedule	Monday through Friday; 8 a.m.	
Expected TAT	2-4 Days (Not reported on Saturday or Sunday)	
Notes	Specimen Stability: Refrigerated (preferred) - 90 Days Ambient - 90 Days Frozen - 90 Days	
CPT Code(s)	82542	
Lab Section	NRLS-Mayo Medical Laboratories	

Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025