Labcorp Oklahoma, Inc. Test Directory

Eastern Equine Encephalitis Antibody, IgG and IgM

Order Name: **EEEP**

Revision Date: 12/05/2019

TEST NAME	METHODOLOGY	LOINC CODE
Eastern Equine Encephalitis Antibody, IgG and IgM		

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.5 mL (0.15 mL)	Serum	Serum Gel Tube	Refrigerated	
Alternate 1	0.5 mL (0.15 mL)	Serum	Red Top	Refrigerated	
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.				
Reference Range	IgG: <1:10 IgM: <1:10 Reference values apply to all age:	s.			
Methodology	Immunofluorescence Assay (IFA)				

GENERAL INFORMATION		
Testing Schedule	Varies	
Expected TAT	1-4 Days	
Notes	Specimen Stability: Refrigerated(preferred) - 14 Days Frozen - 14 Days	
CPT Code(s)	86652 x 2	
Lab Section	NRLS-Mayo Medical Laboratories	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025