

Inhibin A, Tumor Marker

Order Name: **INHA**
Revision Date: 12/09/2019

TEST NAME	METHODOLOGY	LOINC CODE
Inhibin A, Tumor Marker		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.6 mL (0.5 mL)	Serum	Red Top	Refrigerated
Alternate 1	0.6 mL (0.5 mL)	Serum	Serum Gel Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered			
Reference Range	Males: <2.0 pg/mL Females <11 years: <4.7 pg/mL 11-17 years: <97.5 pg/mL Premenopausal: <97.5 pg/mL Postmenopausal: <2.1 pg/mL			
Methodology	Sequential 2-Step Immunoenzymatic Assay			

GENERAL INFORMATION	
Testing Schedule	Monday-Friday 5 a.m.-12 a.m.; Saturday 6 a.m.-6 p.m.
Expected TAT	1-3 Days
Notes	Specimen Stability: Refrigerated(preferred) - 7 Days Frozen - 90 Days
CPT Code(s)	86336
Lab Section	NRLS-Mayo Medical Laboratories