

Organic Acids Screen, Urine

Order Name: **OAU**
Revision Date: 12/09/2019

TEST NAME	METHODOLOGY	LOINC CODE
Organic Acids Screen, Urine		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (4 mL)	Urine, Random	Container, preservative-free	Frozen
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Methodology	Gas Chromatography-Mass Spectrometry (GC-MS)			

GENERAL INFORMATION	
Testing Schedule	Monday-Saturday; 7 a.m.
Expected TAT	3-6 Days
Notes	Specimen Stability: Frozen (preferred) - 416 Days Refrigerated - 14 Days
CPT Code(s)	83919
Lab Section	NRLS-Mayo Medical Laboratories