Labcorp Oklahoma, Inc. Test Directory

Epstein-Barr Virus (EBV) Antibody Profile

Order Name: **SEBV**

Revision Date: 12/13/2019

| TEST NAME | METHODOLOGY | LOINC CODE |
|---|-------------|------------|
| Enctain Barr Virus (EDV) Antibody Brofile | | |

Epstein-Barr Virus (EBV) Antibody Profile

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1.0 mL (0.6 mL) | Serum | Serum Gel Tube | Refrigerated |
| Alternate 1 | 1.0 mL (0.6 mL) | Serum | Red Top | Refrigerated |
| Instructions | Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered. | | | |
| Reference Range | Negative | | | |
| Methodology | Multiplex Flow Immunoassay | | | |

| GENERAL INFORMATION | | |
|---------------------|---|--|
| Testing Schedule | Monday through Friday; Continuous 9 a.m6 p.m., Sunday; 6 a.m. | |
| Expected TAT | 1-3 Days | |
| Notes | Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 14 Days | |
| CPT Code(s) | 86664-EBNA 86665 x 2-VCA, IgG and IgM | |
| Lab Section | NRLS-Mayo Medical Laboratories | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025