

## Epstein-Barr Virus (EBV) Antibody Profile

Order Name: **SEBV**  
Revision Date: 12/13/2019

TEST NAME	METHODOLOGY	LOINC CODE
Epstein-Barr Virus (EBV) Antibody Profile		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>1.0 mL (0.6 mL)</b>	<b>Serum</b>	<b>Serum Gel Tube</b>	<b>Refrigerated</b>
<b>Alternate 1</b>	<b>1.0 mL (0.6 mL)</b>	<b>Serum</b>	<b>Red Top</b>	<b>Refrigerated</b>
<b>Instructions</b>	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
<b>Reference Range</b>	Negative			
<b>Methodology</b>	Multiplex Flow Immunoassay			

GENERAL INFORMATION	
<b>Testing Schedule</b>	Monday through Friday; Continuous 9 a.m.-6 p.m., Sunday; 6 a.m.
<b>Expected TAT</b>	1-3 Days
<b>Notes</b>	Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 14 Days
<b>CPT Code(s)</b>	86664-EBNA  86665 x 2-VCA, IgG and IgM
<b>Lab Section</b>	NRLS-Mayo Medical Laboratories