

Extractable Nuclear Antigen Panel

Order Name: **ENA**  
Revision Date: 04/23/2020

TEST NAME		METHODOLOGY		LOINC CODE
Extractable Nuclear Antigen Panel				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0 mL (1.0 mL)	Serum	Red Top	Refrigerated
Alternate 1	2.0 mL (1.0 mL)	Serum	Serum Gel Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	Multiplex Flow Immunoassay			

GENERAL INFORMATION	
Testing Schedule	Monday-Saturday; 4 p.m.
Expected TAT	1-3 Days
Notes	Panel includes: JO1 AB, Ribonuclear Protein (RNP) AB, SCL70 AB, Smith AB, SSA, SSB AB.
CPT Code(s)	86235x6
Lab Section	NRLS-Mayo Medical Laboratories