

## Extractable Nuclear Antigen Panel

Order Name: **ENA**  
Revision Date: 04/23/2020

TEST NAME	METHODOLOGY	LOINC CODE
Extractable Nuclear Antigen Panel		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>2.0 mL (1.0 mL)</b>	<b>Serum</b>	<b>Red Top</b>	<b>Refrigerated</b>
<b>Alternate 1</b>	<b>2.0 mL (1.0 mL)</b>	<b>Serum</b>	<b>Serum Gel Tube</b>	<b>Refrigerated</b>
<b>Instructions</b>	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
<b>Reference Range</b>	An interpretive report will be provided.			
<b>Methodology</b>	Multiplex Flow Immunoassay			

GENERAL INFORMATION	
<b>Testing Schedule</b>	Monday-Saturday; 4 p.m.
<b>Expected TAT</b>	1-3 Days
<b>Notes</b>	Panel includes: JO1 AB, Ribonuclear Protein (RNP) AB, SCL70 AB, Smith AB, SSA, SSB AB.
<b>CPT Code(s)</b>	86235x6
<b>Lab Section</b>	NRLS-Mayo Medical Laboratories