

Chromium, Urine

Order Name: **Chromi U Rand**

Test Number: 5194545

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Chromium, Urine		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Specimen Type: Plastic urine container, no preservative Specimen Storage: Maintain specimen at room temperature. Specimen Collection: Sampling time is during shift or at the end of shift at the end of the work week for industrial exposure monitoring. For metals with multiphase elimination, the timing is given in relation to the work day exposure (shift) as well as to the work week exposure. Special Instructions: Do not use preservative. Preservatives used for routine analysis may contain mercuric oxide (ie, Stabilur), which interferes with all metal testing. If both urinalysis and metal testing are ordered, please submit a separate urine specimen (containing no a Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days			

GENERAL INFORMATION	
Performing Labcorp Test Code	071530
Notes	Labcorp Test Code: 071530
CPT Code(s)	82495, 82570
Lab Section	Reference Lab