

## Manganese, Urine

Order Name: **Mang U Rand**

Test Number: 5194547

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Manganese, Urine		

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Instructions</b>	<b>Specimen Type:</b> Plastic urine container, no preservative <b>Specimen Storage:</b> Maintain specimen at room temperature. <b>Specimen Collection:</b> Not Available <b>Special Instructions:</b> Do not use preservatives. Preservatives used for routine analysis may contain mercuric oxide (ie, Stabilur), which interferes with all metal testing. If both urinalysis and metal testing are ordered, please submit a separate urine specimen (containing no a <b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days			

### GENERAL INFORMATION

<b>Performing Labcorp Test Code</b>	071597
<b>Notes</b>	Labcorp Test Code: 071597
<b>CPT Code(s)</b>	82570, 83785
<b>Lab Section</b>	Reference Lab