

VMA, Random Urine

Order Name: **VMA U Rand**

Test Number: 5194550

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
VMA, Random Urine		

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	Specimen Type: Plastic urine container Specimen Storage: Vanillylmandelic acid is STABLE at room temperature unpreserved for seven days. Specimen Collection: If specimen is to be stored for longer than seven days before analysis, use random urine transport tube. Special Instructions: Use for patients where 24-hour urine collection is difficult (eg, pediatric patients). Twenty-four hour urine collection is preferable (see Vanillylmandelic Acid (VMA), 24-Hour Urine [004143]). Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days			

GENERAL INFORMATION

Performing Labcorp Test Code	123208
Notes	Labcorp Test Code: 123208
CPT Code(s)	82570, 84585
Lab Section	Reference Lab