

## ToxAssure Flex 14, Urine

Order Name: **TA Flex 14 UR**

Test Number: 5195610

Revision Date: 01/09/2024

TEST NAME	METHODOLOGY	LOINC CODE
ToxAssure Flex 14, Urine	Immunoassay (IA)	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

### GENERAL INFORMATION

Expected TAT	5-8 days
Performing Labcorp Test Code	912336
CPT Code(s)	80347, 80354, 80307, 80348
Lab Section	Reference Lab