

## Specific Gravity Fluid

Order Name: **SP GRAV FL**  
Test Number: 0800900  
Revision Date: 10/31/2019

TEST NAME	METHODOLOGY	LOINC CODE
Specific Gravity Fluid	Optical	2964-5

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>3mL (1 mL)</b>	<b>Body Fluid</b>	<b>EDTA (Lavender Top)</b>	<b>Ambient / Refrigerated</b>
Instructions	<p><b>Note fluid type on Tube and Request.</b> Good for all Body Fluids except for CSF, Semen and Synovial. Invert EDTA tube several times to mix well and deliver to laboratory immediately. <b>Specimen stability is 24hrs refrigerated.</b> Rejection criteria: Specimens received in Syringes, urine cups, monovette tubes, red top tubes, and gel-barrier tubes. Frozen or clotted specimens will be rejected.</p>			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1 Day
CPT Code(s)	84315
Lab Section	Hematology