

Theophylline

Order Name: **THEOPHYLLI**
Test Number: 4005050
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Theophylline	Enzyme Immunoassay	4049-3

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	<p>Peak: Draw 2-4 hours after oral dose. Stability: Ambient 8 hours, Refrigerated 3 months, Frozen 3 months.</p>			

GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.
CPT Code(s)	80198
Lab Section	Chemistry