

## Orange Citrus Allergen IgG

Order Name: **ORANGE IGG**  
Test Number: 3666325  
Revision Date: 12/14/2022

| TEST NAME                  | METHODOLOGY | LOINC CODE |
|----------------------------|-------------|------------|
| Orange Citrus Allergen IgG |             |            |

| SPECIMEN REQUIREMENTS |   |               |                           |                         |
|-----------------------|---|---------------|---------------------------|-------------------------|
| Specimen              | Specimen Volume (min)   | Specimen Type | Specimen Container        | Transport Environment   |
| Preferred             | <b>0.3 mL (0.1 mL)</b>  | <b>Serum</b>  | <b>Clot Activator SST</b> | <b>Room Temperature</b> |
| Instructions          | <b>Specimen Type:</b> Red-top tube or gel-barrier tube<br><b>Specimen Storage:</b> Room temperature<br><b>Specimen Collection:</b> Not Available<br><b>Special Instructions:</b> Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy.<br><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months |               |                           |                         |

| GENERAL INFORMATION          |                                       |
|------------------------------|---------------------------------------|
| Expected TAT                 | 2 - 4 days                            |
| Clinical Use                 | Orange Citrus: <i>Citrus sinensis</i> |
| Performing Labcorp Test Code | 605057                                |
| Notes                        | Labcorp Test Code: 605057             |
| CPT Code(s)                  | 86001                                 |
| Lab Section                  | Reference Lab                         |