

Propoxyphene and Metabolite S/P Quantitative

Order Name: **Propoxyphene S/P**

Test Number: 6900555

Revision Date: 05/17/2021

TEST NAME	METHODOLOGY	LOINC CODE
Propoxyphene and Metabolite S/P Quantitative	Quant Liquid Chromatography-Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL (0.3 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Instructions	<p>Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.</p> <p>Specimen Storage: Room Temperature</p> <p>Specimen Collection: Collect specimen immediately prior to next dose.</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	2-5 days
Performing Labcorp Test Code	084756
Notes	Labcorp Test Code: 084756
CPT Code(s)	80367
Lab Section	Reference Lab