

Streptomycin Antimicrobial Serum Level

Order Name: **Strep Antimicro Srm**
Test Number: 6905151
Revision Date: 01/02/2018

TEST NAME	METHODOLOGY	LOINC CODE
Streptomycin Antimicrobial Serum Level	Quantitative High Performance Liquid Chromatography	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Specimen Type: Red-top tube</p> <p>Specimen Storage: FREEZE</p> <p>Specimen Collection: Serum should be separated from cells within two hours of venipuncture. Send serum in a plastic transport tube before freezing.</p> <p>Special Instructions: Testing referred to National Jewish Health NJHCO#PKSM Was SMH NOTE: Drug dose and frequency, including time of last dose should be indicated on request form.</p>			

GENERAL INFORMATION	
Expected TAT	14-17 days
Performing Labcorp Test Code	815139
Notes	Labcorp Test Code: 815139
CPT Code(s)	80299
Lab Section	Reference Lab