

# Tryptase

Order Name: **TRYPT**  
Revision Date: 10/10/2019

TEST NAME	METHODOLOGY	LOINC CODE
Tryptase		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>0.5 mL (0.2 mL)</b>	<b>Serum</b>	<b>Serum Gel Tube</b>	<b>Frozen</b>
<b>Alternate 1</b>	<b>0.5 mL (0.2 mL)</b>	<b>Serum</b>	<b>Red Top</b>	<b>Frozen</b>
<b>Instructions</b>	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
<b>Reference Range</b>	<11.5 ng/mL			
<b>Methodology</b>	Fluorescence Enzyme Immunoassay (FEIA)			

GENERAL INFORMATION	
<b>Testing Schedule</b>	Monday-Friday; 0900 and 1300
<b>Expected TAT</b>	2 Days
<b>Notes</b>	Stability: Frozen (preferred) - 14 Days Refrigerated - 7 Days
<b>CPT Code(s)</b>	83520
<b>Lab Section</b>	NRLS-Mayo Medical Laboratories