

Dehydroepiandrosterone Sulfate (DHEA-S)

Order Name: **DHES1**
Revision Date: 10/10/2019

TEST NAME		METHODOLOGY		LOINC CODE
Dehydroepiandrosterone Sulfate (DHEA-S)				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.6 mL (0.5 mL)	Serum	Red Top	Refrigerated
Alternate 1	0.6 mL (0.5 mL)	Serum	Serum Gel Tube	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	See interpretive report.			
Methodology	Immunoenzymatic Assay			

GENERAL INFORMATION	
Testing Schedule	Monday-Friday; 0600-0000
Expected TAT	2 Days
Notes	Specimen Stability: Refrigerated (preferred) 14 Days Frozen 30 Days
CPT Code(s)	82627
Lab Section	NRLS-Mayo Medical Laboratories