

Testosterone, Total and Free

Order Name: **TGRP**
Revision Date: 10/23/2019

TEST NAME		METHODOLOGY		LOINC CODE
Testosterone, Total and Free				
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.5 mL (1.0 mL)	Serum	Red Top	Refrigerated
Instructions	Label specimen with Mobilab label or with patient name (first and last), date of birth, date and time of collection, collector initials, and test(s) being ordered.			
Reference Range	An interpretive report will be provided.			
Methodology	FRTST (Free Testosterone): Equilibrium Dialysis TTST (Total Testosterone): Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)			

GENERAL INFORMATION	
Testing Schedule	Monday-Friday 2 p.m.; Saturday and Sunday-Varies
Expected TAT	4 Days
Notes	Specimen Stability: Refrigerated (preferred) - 14 Days Frozen - 60 Days
CPT Code(s)	84402 84403
Lab Section	NRLS-Mayo Medical Laboratories