

## Eosinophil Count, Absolute

Order Name: **EOS CT ABS**

Test Number: 0100050

Revision Date: 10/31/2019

TEST NAME	METHODOLOGY	LOINC CODE
Eosinophil Count, Absolute	Flow cytometry	711-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Alternate 1	1 mL (0.5 mL)	Whole Blood	EDTA (Lavender) Microtainer/Bullet	Refrigerated
Instructions	<p><b>For best results</b>                      Room temperature specimens should be tested within 24hrs, otherwise send Refrigerated. Refrigerated specimens can be tested up to 48hrs.                      Specimens received greater than 48hrs old will be canceled.</p>			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Eosinophilia is found in parasitic diseases, allergic diseases, skin disorders, and certain malignancies.
CPT Code(s)	85004
Lab Section	Hematology