

Flow Cytometry on Peripheral Blood

Order Name: **FLOW PB**
 Test Number: 5582600
 Revision Date: 07/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Flow Cytometry on Peripheral Blood	Flow cytometry	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5mL (3mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Alternate 1	5mL (3mL)	Whole Blood	Lithium Heparin (Dark Green Top / No-Gel)	Room Temperature
Instructions	<p>Deliver to laboratory (flow cytometry) ASAP. DO NOT Centrifuge or Refrigerate. Collect Monday through Friday only! Specimen must be received in the lab section by Friday afternoon or collect sample on Monday. Specimen stability: 48hrs Room Temperature.</p>			

GENERAL INFORMATION

Testing Schedule	Tue-Sat
Expected TAT	2 Days
Clinical Use	Assist in the diagnosis of a lymphoproliferative disease.
CPT Code(s)	Test Dependant
Lab Section	Flow Cytometry