

## Protein Urine Timed

Order Name: **PROT TM U**  
Test Number: 3006150  
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Total Urine Volume		3167-4
Creatinine, Urine, mg/dL	Kinetic Alkaline Picrate (Jaffe)	20624-3
Creatinine Urine Timed Estimated		
Creatinine, Urine, 24 Hour		2162-6
Protein 24 Hour Urine mg/dL	Pyrogallol Red	
Protein Urine Timed Estimated		
Protein 24 Hour Urine mg/24hr		2889-4
Hours Collected - Prompt		30211-7
Volume 24 Urine Prompt		

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, 24-hour	24 hour Urine Container	Refrigerated
Instructions	Timed urine collection. No preservative. Record number of hours and volume in ml on the specimen container. <b>Keep refrigerated.</b> Specimen stability: Ambient 24 hours, Refrigerated 6 days, Frozen 1 month			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful as an aid in evaluating renal function.
CPT Code(s)	84156 (81050 if urine measured by RML)
Lab Section	Chemistry