

## Beta-2, Transferrin

Order Name: **BETA-2 TRA**

Test Number: 3656675

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Beta-2, Transferrin	Immunofixation Electrophoresis	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2mL (1mL)	Fluid	Sterile Screwtop Container	Refrigerated
Instructions	<b>Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens!</b> <b>Notes:</b> 1 mL (Note: This volume Does NOT allow for repeat testing). <b>Specimen Type:</b> Sterile container <b>Specimen Storage:</b> Refrigerated <b>Special Instructions:</b> Testing referred to ARUP TC 50047 <b>Specimen Stability:</b> Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks			

GENERAL INFORMATION	
Expected TAT	4 - 8 days
Clinical Use	This test is to detect CSF in body fluids such as sinus or ear drainage.
Performing Labcorp Test Code	829030
Notes	Labcorp Test Code: 829030
CPT Code(s)	86335
Lab Section	Reference Lab