## Labcorp Oklahoma, Inc. Test Directory

## Bilirubin, Body Fluid

Order Name: **BFTBILI**Revision Date: 01/12/2013

TEST NAME			METHODOLOGY	LOINC CODE
Bilirubin, Body Fluid				
SPECIMEN REQUIR	EMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0 mL (1 mL)	Body Fluid	Serum Gel Tube	See Instructions
Instructions	Place body fluid in serum gel vacutainer tube(s). Spin down and send 2 ml. (minimum volume: 1 mL) of fluid.  Keep refrigerated. Avoid hemolysis. Protect from light.			
	Note:  1. Indicate specimen source on request form.  2. Label specimen with patients name( first and last) or other unique identifier, patient hospital identification number or Social Security Number, date and time of draw, collector initials, and test(s) being ordered.			
Reference Range	An interpretative report will be	provided.		
Methodology	Includes bilirubin indirect			

GENERAL INFORMATION		
Testing Schedule	Monday through Saturday TAT: 2 hours STAT TAT: 1 hour	
CPT Code(s)	82247 - bilirubin, total 82248 - bilirubin, direct	
Lab Section	NRLS-Core Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025