

Lipase, Body Fluid

Order Name: **BFLIP**
Revision Date: 01/12/2013

TEST NAME	METHODOLOGY	LOINC CODE
Lipase, Body Fluid		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	1 mL (minimum volume: 0.5 mL) of body fluid in a mint green-top (lithium heparin) plasma gel tube(s) or a screw-capped, sterile container. Send specimen refrigerated. Maintain sterility and forward promptly. Note: 1. Indicate specimen source on request form. 2. Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of collection, collector initials, and test(s) being ordered.			

GENERAL INFORMATION	
Testing Schedule	Monday through Sunday TAT: 2 hours STAT TAT: 1 hour
CPT Code(s)	83690
Lab Section	NRLS-Core Chemistry