

## Vancomycin, Random

Order Name: **VANCR**  
Revision Date: 01/12/2013

TEST NAME	METHODOLOGY	LOINC CODE
Vancomycin, Random		

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Instructions</b>	Draw blood in a serum gel tube(s). Spin down and send 1 mL (minimum volume: 0.5 mL) of serum refrigerated. Note: Label specimen with patient name (first and last) or other unique identifier, patient hospital identification number or Social Security number, date and time of draw, collector initials, and test(s) being ordered.			

### GENERAL INFORMATION

<b>Testing Schedule</b>	Monday through Sunday TAT: 2 hours STAT TAT: 1 hour
<b>CPT Code(s)</b>	80202
<b>Lab Section</b>	NRLS-Core Chemistry