

## Pediatric Environmental Allergen Panel

Order Name: **A PED EN P**  
 Test Number: 5622600  
 Revision Date: 09/11/2023

TEST NAME	METHODOLOGY	LOINC CODE
Aspergillus fumigatus Allergen	ImmunoCAP	15549-9
Cat Dander Allergen	ImmunoCAP	15609-1
Dog Dander Allergen	ImmunoCAP	15685-1
Alternaria alternata Allergen	ImmunoCAP	15530-9
Cladosporium herbarum Allergen	ImmunoCAP	6075-6
House Dust Mite Allergen	ImmunoCAP	15682-8
Dust Mite (D. farinae) Allergen	ImmunoCAP	15680-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3.6 mL (1.8 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p><b>Recommended for 2 years to 12 years of age.</b>  <b>Specimen Type:</b> Red-top tube or gel-barrier tube,                      Separate Serum from Cells into a screwtop transport container.  <b>Stability Requirements:</b> Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)</p>			

### GENERAL INFORMATION

Expected TAT	3-5 Days
Performing Labcorp Test Code	606470
CPT Code(s)	86003x8
Lab Section	Reference Lab